# **FERNO**°



The PRO 28Z Chair Cot — provides exceptional maneuverability and stability. With its specialized ZFrame™ Advanced Positioning, the PRO 28Z transforms from a cot to a chair so that you — can easily navigate tight hallways, narrow staircases, and elevators — all without transferring the patient. The PRO 28Z is compatible with the FERNO iNLINE® Cot Fastening System for SAE-compliant safety in the ambulance.



Swing-up undercarriage

allows easy roll-in and roll-out unloading



Easily transitions from cot to chair position

with ZFrame Advanced Positioning



**Cross-Compatible** 

When paired with the iNLINE Cot Compatibility Kit, the Pro 28Z is compatible with the SAEcompliant iNLINE Cot Fastening System





TIP**RESIST** ADVANCED STABILITY TECHNOLOGY

### more MANEUVERABILITY



- enables the PRO 28Z to transform from a cot to a chair for ultimate maneuverability in confined spaces
- SAE J3027 Compliant with iNLINE Cot Fastening System for in-ambulance safety
- One-Hand release mechanism quickly and easily switches the Pro 28Z between three locked positions: Chair, intermediate, and cot

Certified SAE J3027 Compliant





<sup>\*</sup> Than the previous model, the 28Z PROFlexx, in the seated position.







# **FEATURES AND BENEFITS**

# ENHANCED PATIENT CARE, COMFORT AND SAFETY

- ZFrame Advanced
   Positioning enables the
   Pro 28Z to transform from a cot to a chair, for ultimate maneuverability in confined spaces.
- Most compact chair cot on the market with a footprint of only 47" x 25" in the chair position.
- One-hand release mechanism quickly and easily switches the Pro 28Z between three locked positions: Chair, intermediate and cot.
- TipRESIST Advanced Stability technology provides 68% more tip resistance in the seated position\*
- Swing-up undercarriage allows easy roll-in loading and roll-out unloading (Oneman loading and unloading is possible when the cot is empty, two operators are always required during patient transports)

- SX Surface Xtenders™
   (Optional accessory) increase the patient surface area by 95%, giving a greater feeling of comfort and stability.
- SAE-compliant Ferno
   X-Restraints are color coded for ease of use

# BUILT FOR LONGEVITY AND SIMPLE OPERATION

- Tested in a controlled environment to 10+ years of service
- All-Metal Frame Construction provides, strength, durability, and performance with heavy patients.
- 2-Year Warranty (Ask your Ferno Representative about our extended service agreements)

#### **CROSS-COMPATIBLE DESIGN**

 iNLINE Compatibility Kit allows you to secure the Pro 28Z to the iNLINE Cot Fastening System, making it SAE compliant. The Pro 28Z is also compatible with the 173 Cot Fastening System

## **SPECIFICATIONS**

Item Number	Base model: 0012826
Length (folded)	81 in / 2063 mm
Length (raised)	76 in / 1925 mm
Width (overall)	25 in / 643 mm
Width (patient surface)	18 in / 448 mm
Height (folded)	8 in / 203 mm
Height	54.5 in / 1384 mm
Depth	47 in / 1202 mm
Weight**	96 lb / 43.5 kg
Weight Capacity	700 lb / 318 kg
Fastening Systems	Ferno iNLINE Cot Fastening System** Ferno 173 Cot Fastening System

## STANDARD FEATURES

- Swing-Down Sidearms
- Ferno X-Restraints
- Two-Piece Locking Mattress
- Backrest Equipment Hook

# **COMPLIANCE**\*\*\*

The **PRO 28Z** has been designed to meet international test standards to increase safety for you and your patient.

- SAE J3027
- BS EN1789
- CAAS-GVS
- KKK-A-1822F
- NFPA 1917

- AMD 004
- IP56
- CE
- EN-1865

# **OPTIONAL ACCESSORIES**



Model 274-S PacRac® Equipment Platform

No. 0818952



SX Surface Xtender™ No. 0822497



**Flip-Out Handles** No. 0822137



Sto-Net<sup>™</sup> Storage Net No. 0822501



Backrest Storage Pouch No. 0374907



**Lower Hard Tray** No. 0822494



NeoMate™ Pediatric Restraint System No. 0314116



Pedi-Mate® + Pediatric Restraint System No. 0314115



KangooFix™ Neonatal Restraint System No. 0374906

\*\*Cot is weighed without mattress and accessories

\*\*\*The Pro 28Z must be used with iNLINE\*\*Cot
Fastening System to be SAE compliant. When used
with the Ferno 173 Cot Fastening System (ratchet-style
fastener) it is not SAE compliant but IS AMD-compliant.

ACCESSORIES NOT PICTURED: iNLINE® Compatibility Kit