

28Z PROFlexx[®] Chair Cot



The 28Z PROFlexx[®] Chair Cot **provides maximum versatility in confining environments**. The 28Z transforms from a cot to a chair so you can easily maneuver through tight hallways, narrow staircases, and small elevators, without transferring the patient. The 28Z's light-weight, superior stability, load capacity, and flexibility make it easy to handle in any situation.

- **Unparalleled Maneuverability** with 48" length in chair position provides the shortest length of any cot on the market for working in confined spaces.
- **Lightest Weight** roll-in style cot at 69 lbs, making it easy to handle in any situation.
- **Swing-Up Undercarriage** allows "one-man[®]" loading and unloading of empty cot.

28Z PROFlexx[®] Chair Cot



Features

- **All-Metal Frame Construction** provides strength, durability and performance with heavy patients
- **Pull-Out Side Lift Handles** provide additional grip points for assistance with heavy patients
- **Ergonomic One-Hand Release** changes cot into one of three locked positions: chair, intermediate, and cot
- **Reflective Labeling** provides night-time visibility for added safety (includes optional personalized graphics)
- **Two 6" Swivel Wheels** and two stationary 6" swivel wheels provide easy rolling and superior steering control
- Available with or without a lower hard tray

Specifications

	MAX	MIN
LENGTH	50 in 2040 mm	75 in 1900 mm
WIDTH	25 in 630 mm	
COT HEIGHT	33 in 840 mm	9 in 220 mm
LOADING HEIGHT	27 in 680 mm	
CHAIR HEIGHT	55 in 1400 mm	
CHAIR DEPTH	48 in 1220 mm	
WEIGHT	82 lbs 37.2 kg	
LOAD LIMIT	700 lbs 318 kg	

Accessories

PATIENT CARE



- **28Z Lead Handle**
No. 0822139



- **Swing Out Lift Handles – 7 in**
No. 0822137



- **Swing Out Lift Handles – 10 in**
No. 0822136

STORAGE ACCESSORIES



- **Pac Rac Equipment Platform**
No. 0818952

Learn more at www.fernoems.com/