

EZ Glide® with PowerTraxx™



70 Weil Way
Wilmington, OH 45177 USA

P 877.733.0911 / 937.382.1451
F 888.388.1349 / 937.382.1191

www.FernoEMS.com

STANDARD FEATURES

- Available in red and blue.
- One inch track to ground clearance
- Extending footrest
- All metal frame construction
- Removable ABS panels
- 5 position lift bar
- One Lithium-ion rechargeable battery with battery meter
- One Battery Charger
- Wireless motor controls
- Electronic brake
- Four inch front swivel wheels
- Extra tall six inch rear wheels

SPECIFICATIONS

Height	37.5 in./953 mm.
Width	20.31 in./516 mm.
Folded Depth	10 in./254 mm.
Weight	56.6 lbs./25.7 kg w/battery
Load Limit	500 lbs./227kgs
Power	28 VDC

Ferno reserves the right to change specifications without notice.

ACCESSORIES

- Additional Batteries
- Additional Charger
- Rear-Lifting Handles
- Headrest
- IV Pole
- O₂ Bottle Holder