

# FERNO Aquaboard Spinal System

The safer way to rescue casualties with suspected spinal injuries from swimming pools



FOR MORE DETAILED  
INFORMATION SEE OVER

## Specifications

### Length

1830 mm

### Width

430 mm

### Weight

8 kg

### Load Limit

182 kg/28st/392lbs

**FERNO** – is a world leading innovator and manufacturer of patient transportation solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **Aquaboard Spinal System** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

call +44 (0) 1274 851 999  
email sales@ferno.co.uk  
www ferno.co.uk

 **FERNO**<sup>®</sup>

# FERNO Aquaboard Spinal System

Safely immobilise a casualty in the water, then support recovery at the poolside

Rescuing a casualty from a swimming pool carries its own distinct risks and challenges – especially when spinal injury is suspected. Rescue equipment needs to be light, buoyant and easily manoeuvrable – and it's vitally important to prevent slipping once the equipment is in place.

**Aquaboard Spinal System** is Ferno's uniquely innovative solution to these issues. Its special design features enable it to float up and under the casualty, so they're safely immobilised whilst still in the water, and supported once they're back on the poolside and in transit to hospital.

Developed in conjunction with the RLSS UK, Aquaboard forms a complete system in itself for swimming pool and leisure centre rescue. It's a simple, highly effective means to minimise the risk of further injury.

- **Easily positioned while in the pool** – smoothly glides in place beneath the casualty
- **Prevents slipping back into the water** – Neoprene Back Support holds the casualty firmly but gently in place
- **Simple to use** – convenient colour-coded adjustable strapping system and concave design to aid positioning
- **No need to reposition the patient after rescue** – they can remain secured on the Aquaboard on the way to hospital and during X-rays
- **Safe and comfortable** – Neoprene Back Support and lightweight adjustable head restraint
- **Easy user familiarisation** – supplied with instructional DVD and patient positioning instructions printed on the neoprene support
- Suitable to fit inside the Model 71, Gazelle, Titan, Titan Tapered and Titan 32 range of Basket stretchers



#### Part Nos:

010776501 Aquaboard Complete Spinal System

#### Spares and Accessories:

038245101 Aquaboard Back Support Pad  
038245200 Aquaboard Head Restraint  
038244600 Aquaboard Body Strap Green Buckle  
038244700 Aquaboard Body Strap Blue Buckle  
038244800 Aquaboard Body Strap Orange Buckle  
038244900 Aquaboard Leg Strap Black Buckle  
038245000 Aquaboard Leg Strap White Buckle  
038244401 Aquaboard Strap System (Set of 5)  
038244501 Aquaboard Strap & Support pad kit (Includes Head & Back Pad, and 5 straps)  
010776400 Wall Mount Bracket for Aquaboard  
081895900 Aquaboard DVD user guide

FERNO - Aquaboard FPC004 10.15

call +44 (0) 1274 851 999  
email sales@ferno.co.uk  
www ferno.co.uk



Ferno (UK) Limited, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.