

Bariatric Solutions



Large patients require added support for their larger mass to provide greater comfort and safe, stable transfers and transport with your cot. **Ferno Bariatric Solutions** allow you to quickly and easily convert standard Ferno cots into a bariatric transport solution, eliminating the need to buy a separate, dedicated stretcher.

- Adapt your cots on-scene for bariatric use
- Solutions are available for the iNX®, POWER X1™, and 35X PROFlexx® Ambulance Cots.



The iNX® with SX Surface Xtenders™ offers a true no-lift bariatric solution.



POWER X1[™] boasts a **700 lb** unassisted lift weight capacity from its lowest position.




35X PROFlexx® has a **1,100 lb load capacity** and 360° turning radius in lowest position.



Bariatric

Solutions



Features

FERNO Bariatric Solutions

 provide the flexibility to safely transfer and transport larger patients with a standard Ferno ambulance cot

Quickly convert to a bariatric cot

 on-scene with the SX Surface Xtender or the LBS Jr. Bariatric Board attachment

The SX Surface Xtender™

 offers 5 locked positions, and widens the cot surface of the iNX and POWER X1

The LBS Jr.™

 provides a narrow profile and folding side frames for greater maneuverability through doorways

When not in use, the LBS Jr. Bariatric board folds

 for easy storage in an ambulance compartment

POWER X1 and 35X PROFlexx cots

 feature 1,100 lb load capacity in the lowest position and true 360° turning radius so you can keep a low center of gravity for moving heavy patients safely

Additional accessories facilitate smooth transfer

 and transport of the patient, including the Manta Transfer Sheet, the Bariatric Tow Package, and Push-Pull Handles

FERNO Bariatric Solutions are Compliant

 with SAE J3027 when combined with corresponding cot fastening system

Solutions & Accessories



SX Surface Extender
 No. FWESX (for iNX)
 No. 1590159 (for POWER X1)



Three-Stage IV Pole
 No. LAIV3SCL (for iNX)
 No. 1590160 (for POWER X1)



 LBS Junior Bariatric Board Includes: LBS Jr., Pillow Wedge, Manta Transfer Sheet, Restraints and Mounting Blocks.
 No. 0822499 (POWER X1)
 No. 0822184 (35X)



LBS Junior Mattress No. 0314050



LBS Junior Carrying Case
 No. 0822179 (mattress only)
 No. 0822180 (mattress & LBS)



Manta Transfer Sheet No. 0382969



IV PoleNo. 0087173 (for LBS Jr.)

Bariatric Tow Package

No. 0822121 (for 35X)



Pillow Wedge No. 0313930



• Straight Push/Pull Handles
No. 0822148
(head end, POWER X1 & 35X)
No. 0822146
(foot end, POWER X1 & 35X)



Angled Push/Pull Handles
 No. 0822147
 (head end POWER X1 & 35X)
 No. 0822112
 (foot end, POWER X1 & 35X)



 *LBS Jr. Adapter Kit for iNX (required for iNX) No. FA911A600

Learn more at www.fernoems.com